

How to measure the following experiments on a GE

1. Preliminary setup

1.1. lock

- for proton set the transmitter frequency around 4.8 ppm *(ppm=4.8)*
- for proton set transmitter power around 300 *(Tx power=300)*
- for a D₂O sample set transmitter power around 150 *(Tx power=150)*
- set the sweep width (sw) to a suitable value
 - to have a relatively large and centered CW signal
 - (e.g. set sw = 1000 Hz
 - adjust the offset
 - set sw = 250 Hz) *(sw = 250 Hz)*
- set the lock gain to see a suitably intense CW signal *(gain = 950)*
- set the lock offset to a value where the CW signal is centered *(offset = -350)*
- set the lock phase to see a nice CW signal
 - starting as a positive peak *(phase = 120)*

1.2. shim

manual shimming: adjust according to the following scheme:

axial z1, z2, z1, z3, z1, (avoid if possible z4 and z5)

radial

back to axial shimming

compu shimming:	<i>initial step</i>	50
	<i>delay</i>	2
	<i>ramp time</i>	2
	<i>AQU</i>	1
	<i>cycles</i>	1

1.3. tuning the probe:

1.3.1. change the relevant cable for the „tune” cable

1.3.2. select the **setup/tune** command from the menu

1.3.3. optimize the signal (as intense and as sharp as possible)

by changing the tune and the match

1.3.4. **quit** form the manual and reconnect the appropriate cable

1.3.5. go back or tune an other channel

2. Set and run a 1D ¹H experiment

2.1. load the *1Puls* pulse sequence (*seqload 1Puls*)

2.2. set F1 channel for proton

<i>(f1 -a)</i>	
<i>nucleus</i>	<i>none</i>
<i>offset</i>	<i>1.6Hz</i>
<i>abs frequency</i>	<i>500.1523816MHz</i>
<i>rel frequency</i>	<i>500.1523800MHz</i>

2.3. set

-the sweep width equ. with 8000 Hz	<i>(sw 8000)</i>
-block size to 4k points	<i>(cb 4k)</i>
-pre acquisition delay (μs)	<i>(de 72)</i>

2.4. adjust also

the gain (use <i>xmt</i> to avoid digitalization overflow)	<i>(gain 3)</i>
the number of scans	<i>(na = 8)</i>
relaxation delay (s)	<i>(rd = 3)</i>
the receiver phase (to avoid zero order phase correction)	<i>(rphase = 30)</i>
the number of dummy goes	<i>(dg = 8)</i>

2.5 run *(zg)*

2.6. process and display (Fourier trans., auto-phase, display) *(ft, aph, dis)*

3. Set for an irradiation (e.g. water suppression)

aim: to place the center of the peak (water peak) on the direct current (DC) signal to be able to use the F1 channel for presaturation.

3.1. Use a narrower spectral width (e.g. $sw = 800$ Hz) to have a higher digital resolution for the same block size (e.g. $cb = 2k$).

3.2. Run the normal 1D experiment.

((zg, ft, aph, dis).

The result is a broad water peak and a sharp direct current (DC) signal.

3.3. Adjust the f1 frequency to have the DC signal on the top (middle) of the H₂O signal
(f1+5h)

memo: if needed return to 3.2.

4. Measure the 90° pulse

aim: to know the length of the time (μs) required at a given power level (e.g. 100% power level) to obtain the bulk magnetization to process completely in the [x, y] plane.

(e.g. seqload 1Puls)

4.1. Measure an approximately 45° (or 60°) 1D proton spectrum.

(pw 5, zg, ft, aph, dis)

4.2. Keep this phase and adjust the power (pw) to obtain now a 360° 1D spec.

(pw 28.8, zg, ft, phase, dis)

4.3. Repeat until no signal intensity left.

4.4. Use the quarter of the final pw level ($1/4 * \text{pw}$) as f190.

5. General settings for 2D experiments

memo: before running the 2D mode of a pulse run the appropriate 1D with a short t_{1dw} (e.g. $t_{1dw} = 10\mu\text{s}$).

before running count the phase cycling steps to set a proper number of scans n_a .

before running guess the value of T_2 . From T_2 calculate $t_{1\text{max}}$.

From $t_{1\text{max}}$ calculate the number of increments when setting up the `acqmode` command.

before running check the total time *(cg et)*

before running check the value of t_{1dw} and make sure that it will reach $t_{1\text{max}}$.

6. Measure a NOESY

Estimate the mixing time.

7. Measure a TOCSY

One has to measure a high power (100% power level) and a low power (e.g. 28% power level) 90° .

The TOCSY mixing time (ms) is:

$t_{mix} * \text{constant} * \text{low power } 90^\circ (\mu s)$

the constant is $4 * 2590 / 90$

rule of thumb for TOCSY transfer in proteins

mixing time (ms) and expected transfer :

$20 < < 30$ from NH to H^α

$30 < < 50$ from NH to H^α and H^β

$70 <$ from NH to H^α , H^β , H^γ and some of the H^δ

To have the effective value of the TOCSY mixing time (ms):

set t_{mix}

run the TOCSY in a 1D mode

type (*cg pv*)

8. Measure a COSY

Run it with more t_1 increments than the appropriate NOESY and/or TOCSY

(e.g. *acqmode 32 0 512 0*)

9. Measure a DQF-COSY

10. Measure a ^1H - ^{15}N HSQC

10.1. calibrate f1 (^1H) as described previously (tune and the 90° as described in 4.)

10.2 calibrate f3 (^{15}N) as below

10.2.1. tune the F3 channel (check the N filter box!)

10.2.2. calibrate the F3 90° pulse length:

(use a suggested value or calibrate with a test sample using a suitable pulse sequence)

10.2.2.1. find a doublet (e.g. an NH proton signal which is a doublet because of the J^{NH} 91Hz)

10.2.2.2. set the power level and the f390 values e.g. (f390=5, and xlev=73)

10.2.2.3. run a 1D and phase it as an anti-phased doublet

10.2.2.4. keep the power level

but alter the f390 values e.g. (f390=43, and xlev=73)

10.2.2.5. keep the values when the doublet vanishes

such a value can be used also for N^{15} decoupling

(e.g. (f390=250, and xlev=53

use the (250*4)ms as value for f3 -g -> f3 -g 1000)

10.3. set (and adjust) the f3 frequency (f3 50.6856770 MHz)

10.4. run the HSQC (seqload glshsqc)

10.5. what are tau=2700, trim1=2000, trim2=1000 values?

11. Measure a ^1H - ^{13}C HSQC

- 11.1. calibrate f1 (^1H) as described previously (tune and the 90° as described in 4.)
- 11.2. calibrate f3 (^{13}C) as below
 - 11.2.1. tune the F3 channel (check the C filter!)
 - 11.2.2. calibrate the F3 90° pulse length:
(use a suggested value or calibrate with a test sample (^{13}C glucose) using a suitable pulse sequence)
 - 11.2.2.1. find the ^1H doublet (the H proton signal which is a doublet because of the ^{13}C : the presence of the J^{CH}).
 - 11.2.2.2. set the power level and the f390 values e.g. (f390= 23.4, and xlev=0)
 - 11.2.2.3. run a 1D and phase it as an anti-phased absorptive doublet
 - 11.2.2.4. keep the f390
but alter the value of the power level e.g. (f390=23.4, and xlev=100)
 - 11.2.2.5. keep the values when the doublet vanishes (zero intensity)
such a value can be used also for C^{13} decoupling
- 11.3. set (and adjust) the f3 frequency (f3 150.92600 MHz)
- 11.4. run the HSQC (seqload gls13ChsqcGrEn_co)
- 11.5. what are tau=1800, tau=900 values?